



Frequently Asked Questions Regarding Life Gear Comfort Zone Blanket Warming System

Answers to questions provided by Casey Jenkins, Life Gear, LLC, casey@life-gear.net

Q) What are the dimensions of your blankets that are included with the warming system, including thickness of each?

A) 40inch x 72inch

Q) What is the construction of the blanket?

A) Waffle weave 100% cotton

Q) What does each blanket weigh?

A) 3 lbs

Q) Is there a warranty on the blankets themselves, separate from the warming system?

A) No. We will replace one if it is received damaged, but that is unlikely it will ever happen. We QC each blanket before putting it in the unit.

Q) How long will the blanket stay warm once it is removed from the box and placed on a patient ?

A) There is really no way to answer this question.. The factors include outside temperature, elevation, weather, movement of patient, etc... Any blanket used will hold heat close to the body, the blanket will reheat as the body expels the heat. This causes the blanket to stay warm for longer periods of time in the areas that are touching the patient or their clothes.

Also, you can use any blanket, towel, or linen of any kind that can handle that heat. They should never use material that has metal or plastic attached. The plates get to 158 degrees, so the blanket will get to the proper temperature of 120 degrees. We have found that very few of our clients use the provided blanket all of the time. It is normally in rotation with a multitude of other blankets.