


# O2-RESQ™

**Pulmonary O2-RESQ™ System offers the clinician the ability to deliver Immediate CPAP Ventilation to the patient.**

**Bilevel Conversion Kit is a cost effective way to convert the O2-RESQ Systems from Immediate CPAP to Bilevel application. The kit includes Exhalation Port, 84" Proximal Pressure Line, Filter, Hose Clips and Optional 24" Tubing with "T".**

The O2-RESQ System is pre-assembled for ease of use and quick set up time. This Single Patient Use O2-RESQ System eliminates the need to keep track of durable equipment and concerns about cross contamination.

Ideal for spontaneously breathing patients suffering from Congestive Heart Failure (CHF) and other Respiratory Distress Conditions.

- *Pre-assembled O2-RESQ System includes:*
  - Flow Generator with Filter
  - 72" Corrugated Anti-Asphyxia Circuit
  - O2-CPAP™ Valve
  - Mask and Head Strap
- *5 Mask styles available: SealMask™, SoftMask™ with Valve, AirMask™, BiTrac ED™ Full Face Mask and BiTrac Shield ED™ Mask*
- *BiTrac ED Full Face Mask features an OmniClip™ with multi-positioning Silicone Forehead Pad*
- *Flow Generator provides 30% O<sub>2</sub> concentration at a Continuous Flow up to 140L/min*
- *O2-CPAP Valves available in: 2.5, 5.0, 7.5, 10.0, 12.5, 15.0 and 20.0cm H<sub>2</sub>O*
- *O2-CPAP Valves provide fixed pressure at any flow rate*
- Single Patient Use
-  Latex Free




*Complete line of pre-assembled CPAP Systems for Immediate/Emergency Ventilation*

# O2-RESQ™

Part No.	Product Description	Size	Qty/Case
<b>Systems</b>			
<b>O2-RESQ System: SealMask and Head Strap, Flow Generator and 72" Corrugated Anti-Asphyxia Circuit with CPAP Valve</b>			
313-7012	SealMask with 10.0cm CPAP Valve	Adult Large	10
313-7011	SealMask with 7.5cm CPAP Valve	Adult Large	10
313-7010	SealMask with 5.0cm CPAP Valve	Adult Large	10
313-7015	SealMask with 10.0cm CPAP Valve	Adult Medium	10
313-7014	SealMask with 7.5cm CPAP Valve	Adult Medium	10
313-7013	SealMask with 5.0cm CPAP Valve	Adult Medium	10
<b>O2-RESQ System: SoftMask with Valve and Head Strap, Flow Generator and 72" Corrugated Anti-Asphyxia Circuit with CPAP Valve</b>			
313-7018	SoftMask with 10.0cm CPAP Valve	Adult Large	10
313-7017	SoftMask with 7.5cm CPAP Valve	Adult Large	10
313-7016	SoftMask with 5.0cm CPAP Valve	Adult Large	10
313-7021	SoftMask with 10.0cm CPAP Valve	Adult Medium	10
313-7020	SoftMask with 7.5cm CPAP Valve	Adult Medium	10
313-7019	SoftMask with 5.0cm CPAP Valve	Adult Medium	10
<b>O2-RESQ System: AirMask and Head Strap, Flow Generator and 72" Corrugated Anti-Asphyxia Circuit with CPAP Valve</b>			
313-7024	AirMask with 10.0cm CPAP Valve	Adult Large	10
313-7023	AirMask with 7.5cm CPAP Valve	Adult Large	10
313-7022	AirMask with 5.0cm CPAP Valve	Adult Large	10
313-7027	AirMask with 10.0cm CPAP Valve	Adult Medium	10
313-7026	AirMask with 7.5cm CPAP Valve	Adult Medium	10
313-7025	AirMask with 5.0cm CPAP Valve	Adult Medium	10
<b>O2-RESQ System: BiTrac ED Full Face Mask with 30/22mm Elbow, OmniClip and Head Strap, Flow Generator and 72" Corrugated Anti-Asphyxia Circuit with CPAP Valve</b>			
313-7035	BiTrac ED Mask with 10.0cm CPAP Valve	Adult Large	10
313-7034	BiTrac ED Mask with 7.5cm CPAP Valve	Adult Large	10
313-7033	BiTrac ED Mask with 5.0cm CPAP Valve	Adult Large	10
313-7038	BiTrac ED Mask with 10.0cm CPAP Valve	Adult Medium	10
313-7037	BiTrac ED Mask with 7.5cm CPAP Valve	Adult Medium	10
313-7036	BiTrac ED Mask with 5.0cm CPAP Valve	Adult Medium	10
<b>O2-RESQ System: BiTrac Shield ED Mask with 30/22mm Elbow and Head Strap, Flow Generator and 72" Corrugated Anti-Asphyxia Circuit with CPAP Valve</b>			
313-7041	BiTrac Shield ED Mask with 10.0cm CPAP Valve	Adult	10
313-7040	BiTrac Shield ED Mask with 7.5cm CPAP Valve	Adult	10
313-7039	BiTrac Shield ED Mask with 5.0cm CPAP Valve	Adult	10
<b>Accessories</b>			
<b>O2-CPAP Valves with Connectors (22mm and 30mm)</b>			
313-7721	20.0cm CPAP Valve (Gold)	313-7717	7.5cm CPAP Valve (Yellow)
313-7720	15.0cm CPAP Valve (Red)	313-7716	5.0cm CPAP Valve (Blue)
313-7719	12.5cm CPAP Valve (Pink)	313-7715	2.5cm CPAP Valve (White)
313-7718	10.0cm CPAP Valve (Green)		
<b>Bilevel Conversion Kit: Exhalation Port, 84" Proximal Pressure Line, Filter, Hose Clips and optional 24" Tubing with "T"</b>			
313-7042			10
<b>Valved "T" Adapter with one arm 22mmF and one arm 22mmM</b>			
313-2061			30
<b>O2-RESQ Filter (30mm M x 30mm F; BFE and VFE 99.9999%)</b>			
313-7043			10
<b>SealMask and Head Strap*</b>		<b>AirMask and Head Strap*</b>	
313-7003	Adult Large	313-7009	Adult Large
313-7002	Adult Medium	313-7008	Adult Medium
313-7001	Adult Small	313-7007	Adult Small
<b>SoftMask with Valve and Head Strap*</b>		<b>BiTrac ED Full Face Mask and Head Strap*</b>	
313-7006	Adult Large	313-7031	Adult Extra Large
313-7005	Adult Medium	313-7030	Adult Large
313-7004	Adult Small	313-7029	Adult Medium
		313-7028	Adult Small
		<b>BiTrac Shield ED Mask and Head Strap*</b>	
		313-7032	Adult

\*All Mask are 10 per case



...bringing change to life™

Pulmodyne  
2055 Executive Drive  
Indianapolis, IN 46241 USA

T: 317.246.5505  
F: 317.246.5501

www.Pulmodyne.com

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

O2-RESQ, O2-CPAP, SealMask, SoftMask, AirMask, BiTrac ED, BiTrac Shield ED, OmniClip, Pulmodyne and bringing change to life are trademarks of Pulmodyne, Inc.